

# Capnography Waveforms – Quick Reference Guide<sup>1</sup>



Normal Capnogram Waveform	Normal Waveform Characteristics	
	<p><b>Phase I:</b> Beginning of exhalation; respiratory baseline containing dead space gas</p> <p><b>Phase II:</b> Rapid upstroke representing exhalation; contains mixture of dead space and alveolar gases</p> <p><b><math>\alpha</math> angle:</b> Normally a 108-degree angle</p> <p><b>Phase III:</b> Alveolar plateau, containing mostly CO<sub>2</sub>-rich alveolar gas.</p> <p><b><math>\beta</math> angle:</b> End of exhalation and beginning of inhalation; normally a 90-degree angle</p> <p><b>EtCO<sub>2</sub>:</b> End of exhaled breath and point of measurement</p> <p><b>Phase IV:</b> Downstroke representing inhalation</p>	<p><b>Normal Range for EtCO<sub>2</sub>:</b></p> <ul style="list-style-type: none"> <li>&gt; 35-45 mmHg</li> <li>&gt; 4.0-5.7 kPa</li> </ul> <p><b>EtCO<sub>2</sub> to PaCO<sub>2</sub> Gradient:</b></p> <ul style="list-style-type: none"> <li>&gt; EtCO<sub>2</sub> is 1-5 mmHg lower than PaCO<sub>2</sub> in patients with normal lung function</li> <li>&gt; Wider gradient indicates greater ventilation/perfusion deficit</li> </ul>
Increasing EtCO <sub>2</sub>	Waveform Characteristics	Nursing Assessment and Considerations
	<ul style="list-style-type: none"> <li>&gt; Increasing amplitude and width, over variable time period, depending on cause</li> <li>&gt; Slowing frequency with decreasing respiratory rate is dependent on cause and patient's physiologic response</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Assess patient for bradypnea/hypercapnia</li> <li>&gt; Assess patient for respiratory failure or over sedation</li> <li>&gt; Consider airway management if needed</li> <li>&gt; Assess for fever or change in temperature from hypothermia to normothermia</li> <li>&gt; Assess for hypermetabolic state</li> <li>&gt; Tourniquet release, sodium bicarbonate, and CO<sub>2</sub> insufflation can cause brief rise in EtCO<sub>2</sub></li> </ul>
Decreasing EtCO <sub>2</sub>	Waveform Characteristics	Nursing Assessment and Considerations
	<ul style="list-style-type: none"> <li>&gt; Decreased amplitude and width</li> <li>&gt; Faster frequency, increased respiratory rate</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Assess patient for tachypnea/hypocapnia.</li> <li>&gt; If tachypneic, assess for underlying causes such as pain, anxiety, or respiratory distress</li> <li>&gt; Assess patient for hypoxemia</li> <li>&gt; Assess patient for decreasing metabolic rate, hypovolemia, or shock</li> <li>&gt; Assess patient for temperature change</li> <li>&gt; Assess patient for pulmonary embolism</li> </ul>
Loss of Waveform	Waveform Characteristics	Nursing Assessment and Considerations
	<ul style="list-style-type: none"> <li>&gt; Loss of capnographic waveform</li> <li>&gt; No breath detected by capnograph</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Assess patient for apnea, complete airway obstruction, or cardiac arrest</li> <li>&gt; If intubated check for ET tube extubation, kinks or blockage, or ventilator disconnection</li> <li>&gt; Confirm cannula or mask is placed on patient correctly and connected to monitor</li> <li>&gt; Ensure patient's airway is open and patent, and patient is breathing</li> <li>&gt; If patient is mouth breathing use cannula with oral prong to capture breaths from mouth</li> <li>&gt; Follow your institution's procedure for airway and breathing support</li> <li>&gt; Check for equipment failure</li> </ul>
Obstructive Airway	Waveform Characteristics	Nursing Assessment and Considerations
	<ul style="list-style-type: none"> <li>&gt; Phase II slopes upward with a blunted <math>\alpha</math> angle instead of a sharp upstroke with strong <math>\alpha</math> angle</li> <li>&gt; Phase III (plateau) is more rounded</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Assess patient for bronchospasm</li> <li>&gt; If intubated, assess ET tube for partial kinking</li> <li>&gt; Assess patient for foreign body in airway</li> <li>&gt; Assess patient for partial airway obstruction</li> <li>&gt; The greater the "shark fin" shape, the greater the severity of the obstructive or reactive airway disease</li> </ul>

<sup>1</sup>Brast, S., Bland, E., Jones-Hooker, C., Long, M., and Green, K. (2016). Capnography for the Radiology and Imaging Nurse: A Primer. *Journal of Radiology Nursing*, Volume 35, Issue 3, 173 - 190. For professional use. See Directions for Use for full prescribing information including indications, contraindications, warnings and precautions. Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

# Capnography Waveforms – Quick Reference Guide<sup>1</sup>

Rebreathing CO <sub>2</sub>	Normal Waveform Characteristics	Nursing Assessment and Considerations
	<ul style="list-style-type: none"> <li>&gt; Waveform shape may be normal but appears to float above baseline</li> <li>&gt; Phase IV (downstroke representing inhalation) does not go back to baseline or zero</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Assess patient for rebreathing of exhaled CO<sub>2</sub></li> <li>&gt; If patient is mechanically ventilated, assess for air trapping or breath stacking, check vent settings, circuit set up, and equipment</li> <li>&gt; If patient is spontaneously breathing, ensure drapes are not covering face</li> <li>&gt; If using oxygen mask, ensure oxygen flow is adequate</li> </ul>
Leak	Waveform Characteristics	Nursing Assessment and Considerations
	<ul style="list-style-type: none"> <li>&gt; First waveform is normal but second and third waveforms (phase III and phase IV) are degraded or distorted</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Assess for mask leak or loss of seal in patients on CPAP, BiPAP or NIV</li> <li>&gt; For intubated patients check ET tube cuff for leak</li> </ul>
Curare Cleft and Secondary (Camel) Hump	Waveform Characteristics	Nursing Assessment and Considerations
	<ul style="list-style-type: none"> <li>&gt; Notching or cleft in first two waveforms during Phase III</li> <li>&gt; Small, secondary waveform during Phase I, between second and third waveform</li> </ul>	<ul style="list-style-type: none"> <li>&gt; Sometimes seen in mechanically ventilated patients due to weak, uncoordinated diaphragmatic movement or ventilator asynchrony under sedation and clearing chemical paralysis</li> <li>&gt; Reassess patient's need for sedation and chemical paralysis</li> <li>&gt; May also be seen with neuromuscular dysfunction</li> </ul>
EtCO <sub>2</sub> Trend Data During Cardiac Arrest		Five-Step Method for Capnography Interpretation
<ul style="list-style-type: none"> <li>&gt; Initially, waveform and EtCO<sub>2</sub> value are low during cardiac arrest but increase with effective chest compressions</li> <li>&gt; The higher the EtCO<sub>2</sub> value the greater the chance of ROSC</li> <li>&gt; EtCO<sub>2</sub> ≥ 20 mmHg during CPR indicates adequate chest compressions</li> <li>&gt; EtCO<sub>2</sub> &lt; 10 mmHg after 20 minutes of high-quality chest compressions indicates poor prognosis for survival</li> </ul>		<ol style="list-style-type: none"> <li><b>1. Is there a waveform? No:</b> Is patient pulseless, apneic, or accidentally extubated? Is patient mouth breathing with a nasal cannula? Is the airway obstructed? Repositioning head may alleviate airway obstruction. <b>Yes:</b> What is the height, width, and frequency (respiratory rate)? Is there a pattern? Is patient hypoventilating or hyperventilating? Is patient in shock or hypovolemic?</li> <li><b>2. What is the shape of the waveform?</b> Do you see a steep rise in Phase II with a plateau? Is there sloping, notching, or a prolonged Phase III? If the plateau is altered, the expiratory phase and alveolar gas exchange are altered. If sloping is seen, consider bronchospasm, kinked artificial airway, or foreign body.</li> <li><b>3. Does the waveform have a steep return to baseline?</b> Phase IV represents the inspiratory phase. If patient is rebreathing CO<sub>2</sub>, assess for air trapping or excess dead space in ventilator circuit.</li> <li><b>4. What is the EtCO<sub>2</sub> trend?</b> Evaluating the trend provides a graphic representation of patient's ventilatory status over time. Downward trending could indicate shock or hypovolemia. Upward trending could indicate increased metabolic demand, hypoventilation, or hyperthermia.</li> <li><b>5. Does your capnographic assessment correlate to your clinical assessment?</b> Are there assessment disparities among EtCO<sub>2</sub> trend, respiratory rate, waveform, and clinical picture? As a standard of practice, multiple subjective and objective assessment criteria are required to confirm ET tube placement.</li> </ol>

<sup>1</sup>Bratt, S., Bland, E., Jones-Hooker, C., Long, M., and Green, K. (2016). Capnography for the Radiology and Imaging Nurse: A Primer. *Journal of Radiology Nursing*, Volume 35, Issue 3, 173 – 190.  
For professional use. See Directions for Use for full prescribing information including indications, contraindications, warnings and precautions. Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.