

03[®]


Regional Oximeter

Masimo Open Connect[®] (MOC-9[®]) Measurement Module


MODULE DIMENSIONS

Typical	Module Dimensions
1.8 in (4.60 cm)	Width
4.0 in (10.2 cm)	Length
1.25 in (3.2 cm)	Thickness
12 ft (3.6 m)	Cable Length

MODULE OUTLINE



INTERFACE CONNECTOR

Connector Description	Power and Communications Connector
Number of Pins	9
Receptacle Part No.	Masimo P/N 27014 or 36694
Connector	MOC-9 [®] Connector 

DC POWER REQUIREMENTS

	Typical	Maximum
Input Voltage	+5V	±0.25V
Line/Load Regulation	-	0.3% rms
Avg. Power	950 mW	1050 mW

COMMUNICATIONS INTERFACE

USB Standard	2.0
USB Data Speed	Full
USB Data Interface	FTDI USB Virtual COM Port (VCP)

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

